

Multiplexed RNA-FISH-guided Laser Capture Microdissection RNA Sequencing Improves Breast Cancer Molecular Subtyping and Prognostic Classification

Evan D. Paul^{1,2}, Barbora Huraiová^{1,2}, Natália Valková^{1,2}, Natália Matyašová^{1,2,3}, Daniela Gábrišová^{1,2}, Soňa Gubová^{1,2}, Helena Ignačková^{1,2}, Tomáš Ondříš^{1,2}, Sílvia Bendíková^{1,2}, Daniel Lovíšek^{1,2}, Michal Gala^{1,2}, Michaela Ristová^{1,2,4}, Iñaki Comino-Méndez⁵, Pavel Morozov⁶, Fresia Pareja⁷, Jakob N. Kather^{8,9,10}, Pavol Čekan^{1,2}

¹ MultiplexDX, s.r.o., Comenius University Science Park, Bratislava, Slovakia.; ² MultiplexDX, Inc., Rockville, MD, USA.; ³ Institute of Clinical Biochemistry and Diagnostics, University Hospital, Faculty of Medicine in Hradec Kralove, Charles University, Hradec Kralove, Czech Republic.; ⁴ Wellcome Centre for Cell Biology, School of Biological Sciences, University of Edinburgh, Edinburgh, Scotland, UK. Electronic address: s1507601@ed.ac.uk; ⁵ Unidad de Gestión Clínica Intercentros de Oncología Médica, Hospitales Universitarios Regional y Virgen de la Victoria. The Biomedical Research Institute of Málaga (IBIMA-CIMES-UMA), Málaga, Spain.; ⁶ Laboratory for RNA Molecular Biology, The Rockefeller University, New York NY, USA.; ⁷ Department of Pathology and Laboratory Medicine, Memorial Sloan Kettering Cancer Center, New York, NY, USA.; ⁸ Else Kroener Fresenius Center for Digital Health, Technical University Dresden, Dresden, Germany.; ⁹ Department of Medicine I, University Hospital Dresden, Dresden, Germany.; ¹⁰ Medical Oncology, National Center for Tumor Diseases (NCT), University Hospital Heidelberg, Heidelberg, Germany.


36P

Background


- Multigene tests provide valuable information about molecular breast cancer (BCa) subtypes (Luminal A, Luminal B, Her2 overexpressed, and basal-like) and prognostic risk groups that differ in terms of prognosis, response to therapy, and clinical outcomes^{1,2}.
- However, multigene tests show only moderate reproducibility at the single-sample level depending on the array platform, tumor composition, gene list, and thresholds³⁻⁵.
- This raises the following questions: *Did I order the right test? Would multiple tests provide better information?*
- Objective** – We aim to assess the level of discordance in multigene tests and determine if combining information from multiple tests improves diagnostic and prognostic performance.

Methods






The mFISHseq (Multiplex8+) assay






Retrospective cohort



Results – Molecular Subtyping



Results – Prognostic risk




Summary and Conclusions

- Like prior reports³⁻⁵, we observed considerable discordance at the single-sample level, which resulted in unstable classification and altered prognosis as the level of discordance increased.
- Approximately 60% of reclassified patients by consensus molecular subtyping received suboptimal treatment.
- By leveraging unbiased transcriptome profiling, the spatially informed mFISHseq (Multiplex8+) assay provides a novel approach to *combine multigene tests, yielding superior diagnostic and prognostic information*. This finding is supported by analysis of the TEAM pathology study⁷.

References: 1. Litton, J. K., et al. *Am Soc Clin Oncol Educ Book* 39, e1–e7 (2019). 2. Kittaneh, M., et al. *Biomark Cancer* 5, 61–70 (2013). 3. Mackay, A. et al. *JNCI: Journal of the National Cancer Institute* 103, 652–673 (2011). 4. Wiegelt, B. et al. *The Lancet Oncology* 11, 339–349 (2010). 5. Barlett, J. M. S. et al. *J Natl Cancer Inst* 108, djw050 (2016). 6. Gendoo, D. M. A. et al. *Bioinformatics* 32, 1097–1099 (2016). 7. Barlett, J.M.S. et al. *npj Breast Cancer* 7, 90 (2021).

Funding and Disclosure

Funding: MultiplexDX and the European Union's Horizon 2020 research and innovation programme under an EIC Accelerator grant (agreement No 946693).
Disclosure: EDP is an employee of MultiplexDX.



Corresponding/Presenting author email address: paul@multiplexdx.com